

D4S234E-DNA segment on chromosome 4 (unique) 234 expressed sequence Gene

Reference: PTXBC001745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D4S234E; D4S234; NEEP21; P21; neuron-specific protein family member 1; brain neuron cytoplasmic protein 1; carboxyterminally EE-tagged neuron-enriched endosomal 21 kDa protein; neuron specific gene family member 1		

For research use only.

GENE DESCRIPTION	DNA segment on chromosome 4 (unique) 234 expressed sequence
GENE ID	27065
ORIGIN SPECIES	Human
ACCESSIONS	BC001745
VECTOR	pDONR223
NCBI SYMBOL	D4S234E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaagttgggaacaatttcagagaagggcaccagcagccgctgctggaggatggcttcgacacattccctgatgacgcccctcgatgcaatcagctgcagttccgccccggataaggtggtcgtgaaactaagaccgagtgatgaacctgaccgcaagaaggga
aagcacgtctccccaaattgctgagttcaccgtcagcatcacggagggtgtcaccgagaggtttaagggtctcgtgttggtctctctcgccctggccttcctcactgcgtctctctggttggtctacaaggtgtacaagtagaccgcgcctgccccgatgggtctcctcaagaacacc
cagtgcatcccaaggcctggagagctactacgcggagcaagactccagtgcccgaggagaattttacacagtcataaaccactacaacctggccaagcagagcatcacgcgtccgtatcgccctggatgtagtctgtcagaagagaagctgtccgagcaggagactgaagcgg
ctgagaagtcagcttag