



Reference: PTXBC001750

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	centrosomal protein of 72 kDa; centrosomal protein 72kDa; centrosomal protein 72		
For research use only.			

GENE DESCRIPTION	centrosomal protein 72kDa
GENE ID	55722
ORIGIN SPECIES	Human
ACCESSIONS	BC001750
VECTOR	pDONR223
NCBI SYMBOL	CEP72
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr