

MEA1-male-enhanced antigen 1 Gene

Reference: PTXBC001754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HYS; MEA; male-enhanced antigen 1; Male-enhanced antigen (H-Y structural gene)		
For research use only.			

GENE DESCRIPTION	male-enhanced antigen 1
GENE ID	4201
ORIGIN SPECIES	Human
ACCESSIONS	BC001754
VECTOR	pDONR223
NCBI SYMBOL	MEA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcctgaaagcactctgtcaggcgcctcgcctggtgacacagtagtctaggaggagacacatgggcctgagcgtatctcccaatcagactgagggaactgggacatcaggccctcagaaggcactgggattggagcagtgaggagcctgaggagagcagga
ggaaacggggctggccagctggctactctaccagccctgaaccaagatcctgaacaaggagggtggaaactggcaccagtgggggatggagatgtgttctgacatccaggatcgaatcaggccctgggcttatttccagaccaccattagagatgaagatgaagat
gaggaggagctacagcgttgaacaaccacagctctattccatggaccagaacatgtagagctggtgaaaaggacaatggctggagtaagcctgcctgcgccaggggttcctgcctgggctcgggagatcggatgccagtgaggagatgtgtacagaaagccctccaagcc
cggcaggcatccctgcctggaagtga