

CLDND1-claudin domain containing 1 Gene

Reference: PTXBC001757

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf4; GENX-3745; claudin domain-containing protein 1; claudin domain containing 1 protein; membrane protein GENX-3745; claudin domain containing 1		

For research use only.

GENE DESCRIPTION	claudin domain containing 1
GENE ID	56650
ORIGIN SPECIES	Human
ACCESSIONS	BC001757
VECTOR	pDONR223
NCBI SYMBOL	CLDND1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataaccgtttgctacagcatttgtaattgcttgtgtgcttagcctcatttccaccatctacatggcagcctccattggcacagactctggtatgaatcgaagtcaggttcaagaaattccagtgattgaataaaagcatctgggatgaattcatttagtgaggcagatgaaaagactt
ataatgatgcacttttgcataaatggcagctgggattgtggagacgggtatcaccatacccaaaacatgcattggtatagcccacagaaaggacagagtcatttgatgtggtcacaaaatgtgtgatttcacactaactgagcagttcatggagaaatttggatcccggaaccac
aatagcgggattgatccttaggacatcttggcgttgccagttcctttacctttgtgagtttaggtttgatgtgctttggggcctttgacggactttgtgcttgattgccgaagcttatccaccattgccacgggcatctccatctcttcaggtctgtgtactgggtcagtaagtgtg
ttatgttgcgggaattgaactactccaccagaaactagagctccctgacaatgatccggtgaatttggatggtcctctgctgctgtgtctctgctccctacagttcatggtctctgctctcttcctctgggctgctcaccaaccggaagagtacacctaataaggcatatcgtgtggc
atga