



Reference: PTXBC001760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AP4; ZNF464; Zfp464; replication initiator 1; 60 kDa origin-specific DNA-binding protein; 60 kDa replication initiation region protein; ATT-binding protein; DHFR oribeta-binding protein RIP60; H_DJ0584D14.12; replication initiation region protein (60kD);		

For research use only.

GENE DESCRIPTION	replication initiator 1
GENE ID	29803
ORIGIN SPECIES	Human
ACCESSIONS	BC001760
VECTOR	pDONR223
NCBI SYMBOL	REPIN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr