

S100B-S100 calcium binding protein B Gene

Reference: PTXBC001766

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NEF; S100; S100-B; S100beta; protein S100-B; S-100 calcium-binding protein, beta chain; S-100 protein subunit beta; S100 calcium-binding protein, beta (neural); S100 calcium binding protein B		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein B
GENE ID	6285
ORIGIN SPECIES	Human
ACCESSIONS	BC001766
VECTOR	pDONR223
NCBI SYMBOL	S100B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgagctggagaaggccatgtgtgccctcatcgacgtttccaccaatattctggaaggaggagacagcacaagctgaagaaatccgaactcaaggagctcatcaacaatgagctttccatttcttagaggaaatcaaagacaggaggtgtggacaaatcatggaaacac
tggacaatgatggagacggcgaatgtgacttcaggaaatcatggcctttgtgcatgtgtactactgcctgccacgagttcttgaacatgagtga