

LSM1-LSM1 homolog, U6 small nuclear RNA associated (*S. cerevisiae*) Gene

Reference: PTXBC001767

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LSM1 homolog, mRNA degradation associated; LSM1-like protein U6 small nuclear RNA associated; LSM1, U6 small nuclear RNA associated; LSM1 mRNA degradation associated; LSM1 homolog, U6 small nuclear RNA associated; U6 snRNA-associated Sm-like protein LSM1;		

For research use only.

GENE DESCRIPTION	LSM1 homolog, U6 small nuclear RNA associated (<i>S. cerevisiae</i>)
GENE ID	27257
ORIGIN SPECIES	Human
ACCESSIONS	BC001767
VECTOR	pDONR223
NCBI SYMBOL	LSM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactatatgcctggcaccgccagcctcatcgaggacattgacaaaaagcacttggtctgcttcgagatggaaggacacttaggcttttaagaagcattgcaatttgcacaaacttagtgctacatcagactgtggagcgtattcatgtgggcaaaaaatacggatattctcagagg
atttttgggtcagaggagaaaaatgtggtcctactaggagaaatagacttggaaaaggagatgacacccctccagcaagtatcattgaagaattctagaagaacaaagggtggaacagcagaccaagctggaagcagagaagtgaagtgcaggccctgaaggaccgaggtct
ttccattctcgagcagatactcttgatgagtactaa