

GTF2F2-general transcription factor IIF, polypeptide 2, 30kDa Gene

Reference: PTXBC001771

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent helicase GTF2F2; BTF4; RAP30; TF2F2; TFIIF; general transcription factor IIF subunit 2; TFIIF-beta; general transcription factor IIF 30 kDa subunit; general transcription factor IIF, polypeptide 2, 30kDa; transcription initiation factor IIF		

For research use only.

GENE DESCRIPTION	general transcription factor IIF, polypeptide 2, 30kDa
GENE ID	2963
ORIGIN SPECIES	Human
ACCESSIONS	BC001771
VECTOR	pDONR223
NCBI SYMBOL	GTF2F2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgagcgcgggaactcgactgaccggcgccaacagacagaggtgtggctagtcgaagttcctaataattgtcacagcaatgggctaaagcttctggaagaggtgaagttgggaaactcgggattgccaagactcaaggaaggactgaggtgtcatttactttgaatgag
gatcttgcataatcatgatatgttggaacacagcttcagtcagtcctagagaacatccattgtcttgcgaagtgttgaggacagacattaacagtatttactgagagctcatcagataagctgtcattggaaggaaatggtgtacaaagagctgaatgccgaccagctgccagtga
aactacatgcgattaaaaagattgcaaatagaagagcttccaacacagtgaggctatcacaacagctggacaagtgtgaacaaccaattacaacctgttgaatcatcaataatcgaatgaaaggaaaaaagaagacggaagcggagctcgagctgataaacacatgt
tttagacatgctattttcagccttggagaacatcaatactataatcttaaggacttggtagacatcacaacaacctgtgggtgtacctgaaggaaatcttaaaagaaattgtgttcgaatgtaaaaggatccacaaaacacatgggagctgaagccagagtagacacactatcaggag
aagaaaagagtactaa