

LOC284912-hypothetical gene supported by BC001801 Gene

Reference: PTXBC001801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical gene supported by BC001801
GENE ID	284912
ORIGIN SPECIES	Human
ACCESSIONS	BC001801
VECTOR	pDONR223
NCBI SYMBOL	LOC284912
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacttctcagccacaaatcacccttcagggttagaccgagagacacctcccttctggtgaaggagagcccgatacagctggtgccggtggccaagaaggaggtggacaaaacctggcgaccgtccctggcaccacctcaccagaggacagctcacacctccacagggg
tcagaagaaggagctgccctgggatgctgtgtatccctacaactgccaggcacatgccccttga