

HIG2-hypoxia-inducible protein 2 Gene

Reference: PTXBC001863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIG2; C7orf68; HIG-2; hypoxia-inducible lipid droplet-associated protein; hypoxia inducible gene 2; hypoxia-inducible gene 2 protein; hypoxia-inducible protein 2; hypoxia inducible lipid droplet associated		

For research use only.

GENE DESCRIPTION	hypoxia-inducible protein 2
GENE ID	29923
ORIGIN SPECIES	Human
ACCESSIONS	BC001863
VECTOR	pDONR223
NCBI SYMBOL	HIG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcatgtgtgaacctctacctgttaggtgtgtgactgacctactctccattctgttagagtgatggagtccttagagggcttactagagagccatgcctgggacctcttgaccaccagaagccaactagccaacacagagccaccaagggccttcagaccatccatccag
aagcatgtga