



Reference: PTXBC001870

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIF-2B; EIF2B; EIF2Bdelta; translation initiation factor eIF-2B subunit delta; eIF-2B GDP-GTP exchange factor subunit delta; eukaryotic translation initiation factor 2B subunit 4 delta; eukaryotic translation initiation factor 2B, subunit 4 delta, 67kDa;		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 2B, subunit 4 delta, 67kDa
GENE ID	8890
ORIGIN SPECIES	Human
ACCESSIONS	BC001870
VECTOR	pDONR223
NCBI SYMBOL	EIF2B4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr