

C20orf24-chromosome 20 open reading frame 24 Gene

Reference: PTXBC001871

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf24; PNAS-11; RIP5; rab5-interacting protein; chromosome 20 open reading frame 24		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 24
GENE ID	55969
ORIGIN SPECIES	Human
ACCESSIONS	BC001871
VECTOR	pDONR223
NCBI SYMBOL	C20orf24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcggcggcggaaggagagccgctcagccgagctggccaacggggccctcaaatgtctcctctggagtaaggctgctggagcgcgcggcctgggaggataaggatgaatttttagatgtgatctactggtccgacagatcattgtgtgctcgtggtgcatttg
gggagttttgccattacgagggttcttgggaatagcaggattctgctgatcaatgcaggagtcctgtacctctacttcagcaattacctacagattgatgaggaagaatatggtggcacgtgggagctcacgaagggaagggtttatgacctctttgacctgttcatgtatgttagctgatagt
ttacaacaggtcatttggatcatcttttacctgccatccattatga