



Reference: PTXBC001872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transcription factor binding to IGHM enhancer 3; RCCP2; RCCX1; TFEA; bHLHe33; transcription factor E3; class E basic helix-loop-helix protein 33; transcription factor E family, member A; transcription factor for IgH enhancer; transcription factor for immu		

For research use only.

GENE DESCRIPTION	immunoglobulin heavy constant mu
GENE ID	3507
ORIGIN SPECIES	Human
ACCESSIONS	BC001872
VECTOR	pDONR223
NCBI SYMBOL	IGHM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr