



Reference: PTXBC001880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CL6; insulin-induced gene 1 protein; INSIG-1 membrane protein; insulin induced gene 1		

For research use only.

GENE DESCRIPTION	insulin induced gene 1
GENE ID	3638
ORIGIN SPECIES	Human
ACCESSIONS	BC001880
VECTOR	pDONR223
NCBI SYMBOL	INSIG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgtccagcattgcacgaccattcttgagctgctctgtcgcacagcgcgagggccgagggcccccgagccagcgcgcggggctgtggcgccaaggttggggagatgalcaactgttccgtgtccgggcccctcctgtcgcgcggccacgggtccccggacgtgaccccgcgccagggggccgcagtgctgctgatgagcgcgcccgagccggagccctaccccaacactggcatcatcgcctgttgcagaggagcctcgtgctcttctcgttgggggtgtgtctagccctgtgtcctaacctcgtcgcgatccagaggaaatgtactctcttccccgaggaggtgatccgaccattttctccgctcgttgggttctcctcctcgtcgtgggacagactgtctgttggcctactgtaccctgtatcagacatcactcggagaaaccccaaaatttagagagaaatggccaggtgtcatgcgctgcgatcagatttttgttggcattaacccgcagctgtctaatggcatttggccaataatgcacgcttgcgaacttaccgacctatttgggctttggcgacattgatctttccacaggtgcctttgggctggggatccacatagctttcgtgctacgtcgtatccgtgtatcgttataatgggtgtctacgtataccccagatttcctctattcttcttgcctgtattatttcttcagagcgtcagcgtgttcgggaacatagcagcagatttgccttgcgaagggcccatagctgattga