

RPL5-ribosomal protein L5 Gene

Reference: PTXBC001882

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA6; MSTP030; PPP1R135; 60S ribosomal protein L5; protein phosphatase 1, regulatory subunit 135; ribosomal protein L5		

For research use only.

GENE DESCRIPTION	ribosomal protein L5
GENE ID	6125
ORIGIN SPECIES	Human
ACCESSIONS	BC001882
VECTOR	pDONR223
NCBI SYMBOL	RPL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtttgttaaaagtgttaagaataaggcctactttaagagataccaagtgaatttagaagacgacgagagggtaaaactgattattatgctcggaaacgcttgatgatacaagataaaataataacacacaccaatacagatgatagttcgtgtgacaaacagagatattgtcagattgcttatgccgtatagaggggatgatgatagctgcgcagcgtatgcacacgaactgccaaaatatggtgtgaagggttgccctgacaaattatgctgcagccaagtggagggtga