

NDUFS5-NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase) Gene

Reference: PTXBC001884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-15k; CI15K; NADH dehydrogenase [ubiquinone] iron-sulfur protein 5; CI-15 kDa; NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase); NADH:ubiquinone oxidoreductase 15 kDa IP subunit; complex I-15 kDa; NADH:ubiquinone oxidore		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase)
GENE ID	4725
ORIGIN SPECIES	Human
ACCESSIONS	BC001884
VECTOR	pDONR223
NCBI SYMBOL	NDUFS5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccttcttgacatccagaaaagggttcggccttaacatagatcgatggtgacaatccagagtggtgaacagccctacaagatggctggtcgatgccatgctttgaaaaaagatggatagaatgtgcacatggaatcggttatactcggcgagagaaagagtgcagatagaatgatgatttcgtagagtggttgccttcgcagaaaacgatgagacgtgcaggtaccatcagggaagcagcgggataaagtgataaaggaaagtaacccctccacctcaccacattggcaagggggagcctcgccctga