

C1orf63-chromosome 1 open reading frame 63 Gene

Reference: PTXBC001892

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0471 protein C1orf63; C1orf63; NPD014; arginine/serine-rich protein 1; arginine and serine rich protein 1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 63
GENE ID	57035
ORIGIN SPECIES	Human
ACCESSIONS	BC001892
VECTOR	pDONR223
NCBI SYMBOL	C1orf63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaactacgtgaacgacatgtggcgggctcggcaggagagaagattcgccctcgacctcggtcggtcggtccagccggtgtgtcgtcgtctaggagccgtcttttcagaagctctcgggccattcccgctcgcagccgttttctcaggagtcggaggagc
aagtccaggtccgttccgaaggcgcaccagcggaagtacaggcgtctactcgggtcactactcgggagccggtcgcgatcccgagccggttaccgagagaggcgtacgggttcaccaggagatactacggtctccttcggttaccggtccggtccgtagcaggtcg
cgctctcggggaaggctgtactcgggaaggcgtacgcatcgccgggacagcgtctactcgggtcttgggtcgcacagtgtacccggaggagcacagcagatggaggacagatccaggacgaggtcgcgggagcagaaccccttctgcttaagtgaagaagatcgaatggagc
tgtagaataagcaaaacacatgcagcgaagctctaggacaacacattgactgacagtagtctcagaactgttcttcagccaaagaacaagccgtgggaatagggtgatcaagtaatgggtgcaaaagcctgaactgtcgaaaaggtaacagaagatggaactcgaatccaa
tggaacacctaccagcaagaagcatagcttttagctctaataattctgtagcaaaagccaatacaaaaatcagctaaagctccacagaagaggcatcttcaagatcaccaaaaatagatcagaaaaaaagtcataatggactgtggataccatctaa