

MPDU1-mannose-P-dolichol utilization defect 1 Gene

Reference: PTXBC001898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDGIF; HBEBP2BPA; Lec35; My008; PP3958; PQLC5; SL15; mannose-P-dolichol utilization defect 1 protein; HBeAg-binding protein 2 binding protein A; suppressor of Lec15 and Lec35 glycosylation mutation homolog; mannose-P-dolichol utilization defect 1		

For research use only.

GENE DESCRIPTION	mannose-P-dolichol utilization defect 1
GENE ID	9526
ORIGIN SPECIES	Human
ACCESSIONS	BC001898
VECTOR	pDONR223
NCBI SYMBOL	MPDU1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggagggcggacggaccgttaaacggctgctggtccgattctttacgtgagaaatgctacgaccaacttttcgttcagtgaggactgtctcagctccctgctcaagattctcctcagcaaaaggcctggggctgggcattgtggctggtcacttctagtaaagctgccccagg
tgtttaaatcctgggagccaagagtgctgaagggttgagtcaccagctctgtaatgctggagctagtggaactgactgggacattggtctacagcactcaaaactccattcagctctgggtgaagccttatctctgagctccagacgatcaccatctgctcctggtcatgcactacag
aggacagactgtgaaaggtgtcgtcttctcgtctgctacggcctggtcctgctggtgcttctcactctgacgcccttgactgtagtcaccctgctccaggcctccaatgtgctgctgtgtgtgggtgggaggtctccaggcagccaccaactaccacaaggccacagggccagct
ctcagccatcacagttctctgctgtttgggggtccctggccgaattctcacttcattcaggaaacggagatccccgatggctgggaccttgtgtgtctctctctgcaacggcctcattcgccaccagctgctcttctactggaatgcaagcctccccacaagcaaaaaggcg
cagtag