



Reference: PTXBC001918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DSAEC; 3-ketoacyl-CoA thiolase, mitochondrial; acetyl-Coenzyme A acyltransferase 2; beta ketothiolase; mitochondrial 3-oxoacyl-CoA thiolase; mitochondrial 3-oxoacyl-Coenzyme A thiolase; acetyl-CoA acyltransferase 2		

For research use only.

GENE DESCRIPTION	acetyl-Coenzyme A acyltransferase 2
GENE ID	10449
ORIGIN SPECIES	Human
ACCESSIONS	BC001918
VECTOR	pDONR223
NCBI SYMBOL	ACAA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr