



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

BAG1-BCL2-associated athanogene Gene

Reference: PTXBC001936

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HAP; RAP46; BAG family molecular chaperone regulator 1; BCL2-associated athanogene; Bcl-2 associating athanogene-1 protein; Bcl-2-binding protein; glucocorticoid receptor-associated protein RAP46; receptor-associated protein, 46-KD; BCL2 associated athanoge

For research use only.

GENE DESCRIPTION	BCL2-associated athanogene
GENE ID	573
ORIGIN SPECIES	Human
ACCESSIONS	BC001936
VECTOR	pDONR223
NCBI SYMBOL	BAG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagaagaaaccgccgccctcgaccggagcgaggatgtaccgggagcgaggagttgacctgagtgaggaaagcgacctggaagtgaaggcgaccgagtgaggaggcgaccaggcggaagagatgaatcgagccaggagtgaccgggacgaggatc
gaccgggagcgaggaggtgaccaggaggaaatgcccgcagctggggctaccgtgactgcaccacagcaatgagaagcagcactctatgttacctccagcaggcgacagtgaaacagttgtccaaagcctggcccgagttgttgaaagggtcataggggttccacagctcttt
cagaactcatattaaagggaatactctgaaggaaatggaacacacgctgtgcagacttggaaatacaagatgtgtccggctcatgttaattgggaaaaagacactccacgaagagggttgaactaaagaagtgtgaacacatttggagaagctctgtagagaatagctgaccacgtgga
agagttgataaagagcttactggaatccaggggttttctgcccaaggatgttccaagctgaagctgacgtctgacactgtatggagcagataaaagccacaatagacagctttagaagcacttggaggagattgacacatgatctgccagaaattcaagaagcagattgaaaggaa
agcctctgtaaaaagagctcagcactccttagcggatgtgacacatgfcagcagacaacatctccagcagctgaacactgggctcagctctacaacatttgcctctcccgagtgta

More info: Please contact us for any question: contact@proteogenix.fr