

BDH2-3-hydroxybutyrate dehydrogenase, type 2 Gene

Reference: PTXBC001953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFA6R; PRO20933; SDR15C1; UCPA-OR; UNQ6308; 3-hydroxybutyrate dehydrogenase type 2; R-beta-hydroxybutyrate dehydrogenase; dehydrogenase/reductase SDR family member 6; oxidoreductase UCPA; short chain dehydrogenase/reductase family 15C member 1; 3-hydroxyb		

For research use only.

GENE DESCRIPTION	3-hydroxybutyrate dehydrogenase, type 2
GENE ID	56898
ORIGIN SPECIES	Human
ACCESSIONS	BC001953
VECTOR	pDONR223
NCBI SYMBOL	BDH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcgaacttgatgggaaagtcacatcctgacggccgtgctcaggggattggccaagcagctgccttagcttttgaagagaagggtgccaagtcatagccacagacattaatgagtcacaaactcaggaaactgaaaagtaccgggtattcaaaactgtgtccttgatgcacaaag
aagaacaaattgatcagtttgccagtgaaagttgagagacttgatgttctctttaatgttgctgttttgcacatcaggaactgtcctggattgtgaggagaaagactgggacttctcagatgaatcctcaatgtgcgcagcatgtacctgatcaaggcattcttctctaaatgcttgctcagaat
ctggcaatattatcaacatgtctctgtggttcagcgtcaaggagttgtgaacagatgtgtgtacagcacaaccaaggcagccgtgattggcctcacaaatctgtgctgcagatttccagcagggcatcaggtgcaactgtgtgtgcccaggaacagttgatacgccatctctaca
agaaagaatacagccagaggaaatcctgaagaggcacggaatgattctgaagagacaaaagacgggaagattcgcaactgcagaagaaatagccatgctctgctgtatttggcttctgatgaatctgcttatgtaactgtaacccgtgcatcattgatggaggctggagctgtga