



Reference: PTXBC001962

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EXDL2; C14orf114; exonuclease 3'-5' domain-containing protein 2; exonuclease 3'-5' domain-like 2; exonuclease 3'-5' domain-like-containing protein 2; exonuclease 3'-5' domain containing 2		

For research use only.

GENE DESCRIPTION	exonuclease 3'-5' domain-like 2
GENE ID	55218
ORIGIN SPECIES	Human
ACCESSIONS	BC001962
VECTOR	pDONR223
NCBI SYMBOL	EXDL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]