

NRSN2-neurensin 2 Gene

Reference: PTXBC001963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf98; dJ1103G7.6; neurensin-2; neurensin 2		
For research use only.			

GENE DESCRIPTION	neurensin 2
GENE ID	80023
ORIGIN SPECIES	Human
ACCESSIONS	BC001963
VECTOR	pDONR223
NCBI SYMBOL	NRSN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccgagctgcaatcgttctctgagctgagccgggccccagcgtggagatggcaagtggatgggtccgctcctacctgcacctcttctatgaggactgtgcaggcactgctctcagcagcagacctgaggaccccggtctgtgcccccgccggccctggccctcact
gtgttggaagatcagcctgtctcgggacccctgcttctgctgctgggtgtggtggtcctgacctggctatgagtgcccccaagctggaggcctcgggtgaggtgagttctgtgttgatcagcgggcagccgactacaaccaggccctgggacactgtcgttggcaggca
cagcgtctgtgtggcagctgagttctgctgccatctgctcttctgggccaatgataggctggctgagccaggacaccaaggcagagcccttgacccccgaagccgacagccagctggaggtcttctggggatgagccagagcagcagttgtcacccatttccgcaatgccagtggc
cagtcattgttctgccacccgagcccttgggcaatcttctgtgcagactatccagcccaaggaggactctga