



Reference: PTXBC001976

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRP21; PRPF21; SAP114; SF3A120; splicing factor 3A subunit 1; SAP 114; pre-mRNA processing 21; pre-mRNA splicing factor SF3a (120 kDa subunit); spliceosome-associated protein 114; splicing factor 3 subunit 1; splicing factor 3a, subunit 1, 120kD; splicing		

For research use only.

GENE DESCRIPTION	splicing factor 3a, subunit 1, 120kDa
GENE ID	10291
ORIGIN SPECIES	Human
ACCESSIONS	BC001976
VECTOR	pDONR223
NCBI SYMBOL	SF3A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]