

TPD52L1-tumor protein D52-like 1 Gene

Reference: PTXBC002375

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D53; tumor protein D53; tumor protein D52-like 1		

For research use only.

GENE DESCRIPTION	tumor protein D52-like 1
GENE ID	7164
ORIGIN SPECIES	Human
ACCESSIONS	BC002375
VECTOR	pDONR223
NCBI SYMBOL	TPD52L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcgcaggcacaaggtttgttgagactgaaccgttgcaagggaacagacgaagatgcagtagccagtgcctgactctctagcatgctctctgaggaggaaaaggaagagttaaagcagagtagttcagctagaagacgaaattacaactacgacaagttttgcagcgaaa
gaaaggcatctagttgagataaaacaaaactcggcatgaacctgatgaatgaattaaacagaacttcagcaaaagctggcatgacatgcagactaccactgcctacaagaaaacacatgaaacccctgagtcacgcaggggcaaaaggcaactgcagcttcagcaacgttggaacggc
catcagcaagaagttcggagacatgagttactccattcgccattccataaagtatgcctgctatgagacgaaagttag