



Antibodies and Proteins

## Espace Européen de l'Entreprise

67300 Schiltigheim, France

03 90 20 54 70 - Fax +33 (0)3

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

**ATP5D-ATP synthase, H<sup>+</sup> transporting, mitochondrial F1 complex, delta subunit Gene**

Reference: PTXBC002389

|                |                                                                                                                                                                                                                                                              |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ATP synthase subunit delta, mitochondrial; F-ATPase delta subunit; mitochondrial ATP synthase complex delta-subunit precursor; mitochondrial ATP synthase, delta subunit; ATP synthase, H <sup>+</sup> transporting, mitochondrial F1 complex, delta subunit |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ATP synthase, H <sup>+</sup> transporting, mitochondrial F1 complex, delta subunit                |
| <b>GENE ID</b>          | 513                                                                                               |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC002389                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ATP5D                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)