



ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

LGALS3BP-lectin, galactoside-binding, soluble, 3 binding protein Gene

Reference: PTXBC002403

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	90K; BTBD17B; CyCAP; M2BP; MAC-2-BP; TANGO10B; gp90; galectin-3-binding protein; L3 antigen; MAC2BP; Mac-2-binding protein; basement membrane autoantigen p105; lectin galactoside-binding soluble 3-binding protein; mac-2 BP; serum protein 90K; transport an		

For research use only.

GENE DESCRIPTION	lectin, galactoside-binding, soluble, 3 binding protein
GENE ID	3959
ORIGIN SPECIES	Human
ACCESSIONS	BC002403
VECTOR	pDONR223
NCBI SYMBOL	LGALS3BP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr