

PEA15-phosphoprotein enriched in astrocytes 15 Gene

Reference: PTXBC002426

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMAT1; HUMMAT1H; MAT1; MAT1H; PEA-15; PED; astrocytic phosphoprotein PEA-15; 15 kDa phosphoprotein enriched in astrocytes; homolog of mouse MAT-1 oncogene; mammary transforming gene 1, mouse, homolog of; phosphoprotein enriched in diabetes; phosphoprotein		

For research use only.

GENE DESCRIPTION	phosphoprotein enriched in astrocytes 15
GENE ID	8682
ORIGIN SPECIES	Human
ACCESSIONS	BC002426
VECTOR	pDONR223
NCBI SYMBOL	PEA15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgagtacgggaccctctgcaagacctgaccaacaacatcaccttgaagatctagaacagctcaagtcggcctgcaaggaagacatccccagcgaagagtgaggagatcactactggcagtgccctggtttagcttctggagagccacaacaagctggacaaagacaacct
ctctacattgagcacatcttggagatctcccgctcctgacctactcactatggtggtgactacagaaccggtgctggaagatctctgaggagatgagctggacaccaagctaacccgtatccccagtgccaagaagtacaaagacattatccggcagccctctgaggaaagatcat
caaatggctccccaccgaagaaggcctga