



Reference: PTXBC002431

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B4Gal-T2; B4Gal-T3; beta4Gal-T2; beta-1,4-galactosyltransferase 2; UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 2; UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase 2; UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 2; UDP-galactose:		

For research use only.

GENE DESCRIPTION	UDP-Gal:betaGlcNAc beta 1,4- galactosyltransferase, polypeptide 2
GENE ID	8704
ORIGIN SPECIES	Human
ACCESSIONS	BC002431
VECTOR	pDONR223
NCBI SYMBOL	B4GALT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]