

C14orf1-chromosome 14 open reading frame 1 Gene

Reference: PTXBC002444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERG28; NET51; chromosome 14 open reading frame 1		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 1
GENE ID	11161
ORIGIN SPECIES	Human
ACCESSIONS	BC002444
VECTOR	pDONR223
NCBI SYMBOL	C14orf1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccgtttcctgaatgtgttaagaagttggctggttatgggttcacatagccatggggaacagctgcagagcttcgagaccacactttctctatgaaaagctctacactggcaagccaaacctgtgaatggcctccaagctcggacctttgggatctggacgtgctctcatcagtg
attcgtgcctctgtgacattgacattcacaacaagacgctctatcacatcacactctggaccttctctctgcccctgggcatcttctctctgagttgtttgtctatggaactgcagctcccacgattggcgtctcggcaccctgatgggtggcaagtttctccatctgggtatgctggtcgggctc
cgggtatctagaagtagaaccagtatccagacagaagaagagaaactga