

EFHD1-EF-hand domain family, member D1 Gene

Reference: PTXBC002449

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MST133; MSTP133; PP3051; SWS2; EF-hand domain-containing protein D1; EF-hand domain-containing protein 1; swiprosin-2; EF-hand domain family member D1		

For research use only.

GENE DESCRIPTION	EF-hand domain family, member D1
GENE ID	80303
ORIGIN SPECIES	Human
ACCESSIONS	BC002449
VECTOR	pDONR223
NCBI SYMBOL	EFHD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagtgaggagctggcgtgcaagctggagcggcgtggcggcgaggagggcggagagagtgccccagctggctccctcggcgccccagccccggagccaagccccgagccgagcctccgcccgtgccccacggccagcggcgacggcgagctgagcgc
ccagctgagccggcgctggacatcaacgagggcgtggcggccccggcgtgcaagggtttcaacccctacacggagttccccggagttcagccggcctcatcaaggacctggagagcatgttcaactgtatgacgtggcgggatggcttcacgacctgagggagctgaa
gctgatgatgagaagctggggccccccagaccacctggcctgaagagcatgatcaaggaggtggatgaggacttcgatggcaagctcagcttccgggagttcctgctcatttccacaaggccggcgaggaggagctgaggaggacagtggtgctgctggcctggcgaagc
tttctgagatcgtatggccctggagggtgtcaagggtgccagaactctttgaagccaaggtccaagcctgtcctcggcagtaagttgaagcagagttgaaagctgagcaagatgagcggagcgggaggaggaggagggcggctccgacgagcagcctccagaactcaa
ggccaacttcaatacatag