

DERL1-Der1-like domain family, member 1 Gene

Reference: PTXBC002457

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DER-1; DER1; derlin-1; DERtrin-1; Der1-like domain family, member 1; degradation in endoplasmic reticulum protein 1; derlin 1		

For research use only.

GENE DESCRIPTION	Der1-like domain family, member 1
GENE ID	79139
ORIGIN SPECIES	Human
ACCESSIONS	BC002457
VECTOR	pDONR223
NCBI SYMBOL	DERL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggacatcggagactggtcaggagatccggcgatcacgcgtattggtcgcgcaccgtcgcgtgcccttggtcggcaaacctcgcctcatcagcccgccctctctgcccgaagccttcttatcgtttcagattggaggccaatcactgccacctttatttc
cctgtgggtccaggaaactggattctttatttggcaatttatattcttatatcagtaattctacgcgacttgaacaggagcttttgatggaggccagcagactatttatcatgctcctcttaactggattgcatcgtgattactggcttagcaatggatagcattgctgatgattcctctgatcatg
tcagtactttatgtctgggcccagctgaacagagacatgatgtatcatttgggttgaacacgatttaaggcctgctatttaccctgggttatccttggaattcaactatcatcggaggctcggtaacatgagcttattggaatctggttggacatctttatttttctaatgttcagataccaat
ggacttggaggagaagaattttctatccacacctcagttttgtaccgtggtcgcctagtaggagaggagatcaggatttgggtgtgccccctgctagcatgaggcgagctgctgatcagaatggcggaggcgaggagacacaactggggccaggccttcgacttggagaccagtga
a