



Reference: PTXBC002469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MENT; protein MENT; methylated in normal thymocytes protein; chromosome 1 open reading frame 56		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 56
GENE ID	54964
ORIGIN SPECIES	Human
ACCESSIONS	BC002469
VECTOR	pDONR223
NCBI SYMBOL	C1orf56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr