

CPLX1-complexin 1 Gene

Reference: PTXBC002471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPX-I; CPX1; complexin-1; CPX I; complexin I; synaphin 2; complexin 1		
For research use only.			

GENE DESCRIPTION	complexin 1
GENE ID	10815
ORIGIN SPECIES	Human
ACCESSIONS	BC002471
VECTOR	pDONR223
NCBI SYMBOL	CPLX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttggatgaagcaggcttagggggccaccaaggacatggggaagatgctgggggtgacgaggagaaggaccagacgcccaagaaggaggaggagcgggcaggggcgctgcgcaggcgaggaggagcgcaaggccaagtacgccaagatggagg
cggagcgcgagggccgtgcgcaggggcatccgagacaagtacggcatcaagaagaaggaggagcgcgaggccgaggccagggcccatggaggccaactccgagggggagcttgacgcggccaagaaggccatcccgccgggctgcggggacgaggtggaggaggagg
acgagagcatcctggacaccgtcatcaagtacctgccgggcccgtgcaggacatgctcaagaagtag