

PFN1-profilin 1 Gene

Reference: PTXBC002475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALS18; profilin-1; epididymis tissue protein Li 184a; profilin I; profilin 1		
For research use only.			

GENE DESCRIPTION	profilin 1
GENE ID	5216
ORIGIN SPECIES	Human
ACCESSIONS	BC002475
VECTOR	pDONR223
NCBI SYMBOL	PFN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgggtggaacgcctacatgcacaacctcatggcgacgggacctgtcaggacgcggccatcgtgggtctacaaggactgcctccgtctggccgcgcgtcccggaacgttcgtcaacatcacgccagctgaggtgggtgtcctggttgcaagaccgggtcaagtttt
acgtgaatgggtgacacttggggccagaaatgttcggtgatccgggactcactgctgcaggatggggaattagcatggatcttctgaccaagagcaccggtggggcccccaccttcaatgtcactgtcaccaagactgacaagacgtagctcctgctgatgggcaagaaggtgtcc
acgggtggtttgatcaacaagaatgttatgaatggcctcccaccttcggcgttcccagttactga