



Reference: PTXBC002492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRM1; tRNA (guanine(26)-N(2))-dimethyltransferase; N(2),N(2)-dimethylguanosine tRNA methyltransferase; TRM1 tRNA methyltransferase 1 homolog; tRNA 2,2- dimethylguanosine-26 methyltransferase; tRNA methyltransferase 1 homolog; tRNA(guanine-26.N(2)-N(2)) met		

For research use only.

GENE DESCRIPTION	TRM1 tRNA methyltransferase 1 homolog (S. cerevisiae)
GENE ID	55621
ORIGIN SPECIES	Human
ACCESSIONS	BC002492
VECTOR	pDONR223
NCBI SYMBOL	TRMT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr