

TMEM14C-transmembrane protein 14C Gene

Reference: PTXBC002496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf53; HSPC194; MSTP073; NET26; bA421M1.6; transmembrane protein 14C		
For research use only.			

GENE DESCRIPTION	transmembrane protein 14C
GENE ID	51522
ORIGIN SPECIES	Human
ACCESSIONS	BC002496
VECTOR	pDONR223
NCBI SYMBOL	TMEM14C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggacactggctcagtagtgcctttgcatgggttggcttggctacgcagcactgggtgcttctgggtggatcattggctatgtaaaagcaggcagcgtgccgtccctggctgcagggtgctctttggcagctagccggcctgggtgcttaccagctgtctcaggatccaaggaacgt
ttgggttttctagctacatctggtaccttggctggcattatgggaatgaggttctaccactctggaaaattcatgcctgcaggtttaattgcaggtgccagtttgcctgatggctgccaaagtggagtttagtatgttcaacagacccattag