



Reference: PTXBC002523

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCP; HCPH; HPTP1C; PTP-1C; SH-PTP1; SHP-1L; SHP1; tyrosine-protein phosphatase non-receptor type 6; hematopoietic cell phosphatase; hematopoietic cell protein-tyrosine phosphatase; protein-tyrosine phosphatase 1C; protein-tyrosine phosphatase SHP-1; prote		

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 6
GENE ID	5777
ORIGIN SPECIES	Human
ACCESSIONS	BC002523
VECTOR	pDONR223
NCBI SYMBOL	PTPN6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]