



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KIAA0494-KIAA0494 Gene

Reference: PTXBC002525

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

KIAA0494; EF-hand calcium-binding domain-containing protein 14; EF-hand calcium binding domain 14

For research use only.

GENE DESCRIPTION	KIAA0494
GENE ID	9813
ORIGIN SPECIES	Human
ACCESSIONS	BC002525
VECTOR	pDONR223
NCBI SYMBOL	KIAA0494
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr