

HSPA4-heat shock 70kDa protein 4 Gene

Reference: PTXBC002526

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APG-2; HEL-S-5a; HS24/P52; HSPH2; hsp70; hsp70RY; heat shock 70 kDa protein 4; epididymis secretory sperm binding protein Li 5a; heat shock 70-related protein APG-2; heat shock 70kD protein 4; heat shock 70kDa protein 4; heat shock protein, 110 kDa; hsp70		

For research use only.

GENE DESCRIPTION	heat shock 70kDa protein 4
GENE ID	3308
ORIGIN SPECIES	Human
ACCESSIONS	BC002526
VECTOR	pDONR223
NCBI SYMBOL	HSPA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggtggtggcatagacctgggtccagagctgctacgtcgtgtggcccgccggcgccgcatcgagactatcgtatgagtatagcgaccgctgcacgcccgtgtcattcttttggctcctaagaatcgttaattggagcagcagctaaaagccaggttaatttaataatgcaaag
aacacagtcgaaggatttaaaagattccatggccgagcattctctgatccattttgtggagcgagaaaaatctaacttgcataatgatattgtgcagttgctacaggattaacaggtataaagggtgacatatatggaggaagagcagaatggaccagtggaatggacaaggagacaacccagg
ccccaggctgctgagcagggtacagacacagctgtgccttcgattcagacaagaagctctcgaatggacattgattga