

ZNF593-zinc finger protein 593 Gene

Reference: PTXBC002580

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZT86; zinc finger protein 593; zinc finger protein T86		
For research use only.			

GENE DESCRIPTION	zinc finger protein 593
GENE ID	51042
ORIGIN SPECIES	Human
ACCESSIONS	BC002580
VECTOR	pDONR223
NCBI SYMBOL	ZNF593
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggcgaagcggcggcgccgactggatgagattcaccgagctgcggcctcaggatccgcacgacccagcccgacccaaagccgagttcgaccccgacctgccagggggcggtctgcaccgctgtctggcctgcgcgaggtactcatcgattccaccaacctgaa
gacccacttcgatccaaagaccacaagaaaaggctgaagcagctgagcgtcgagccctacagtcagggaagagcgggagagggcagcgggtatgggatcctatgtgccccccaggcgctggcagtgccacggaagtgtccactgaggtccctgagatggatacctctaccta