



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD99-CD99 molecule Gene

Reference: PTXBC002584

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD99 molecule; CD99 antigen; HBA71; MIC2; MIC2X; MIC2Y; MSK5X; E2 antigen; MIC2 (monoclonal antibody 12E7); T-cell surface glycoprotein E2; antigen identified by monoclonal 12E7, Y homolog; antigen identified by monoclonal antibodies 12E7, F21 and O13; ce		

For research use only.

GENE DESCRIPTION	CD99 molecule
GENE ID	4267
ORIGIN SPECIES	Human
ACCESSIONS	BC002584
VECTOR	pDONR223
NCBI SYMBOL	CD99
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcgcggggcgtgcgtgcgcgtcgtccttcgcgcctgcgtgggttctcgtcgtccgcgcccgatgggtgtttcgaattatctgatgccttctgcacatgaaaaacaagaacccactgcgaatccccagaacccagctcgtgggagctgcttgaacttaggagatcgtgtttgatg
gagaaaatgacgaccacgaccaccgaaccaccaaacccatgatgcacaatccaaaccccaaccaccttagttccctcgtgtagcttttcagatcgtcgaacttgcggatggcgtttcagcttggtgagaaaggaaaggagcagtcgatgggtggagcgacgccccagggaaagggggaagag
gccgacgccccggcgtgatccccgggattgtggggcgtcgtcgtgcgccgtggtggagccgactcttcagatccttactcaccagaaaagaagctatgcttcaaaagaaatcgacacaaggggaggtggacatggagagccaccgggaatgccaacgcgagccagctgttcagc

More info: Please contact us for any question: contact@proteogenix.fr