

MTCP1-mature T-cell proliferation 1 Gene

Reference: PTXBC002600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P13MTCP1; p8MTCP1; protein p13 MTCP-1; MTCP-1 type B1; mature T-cell proliferation 1 isoform p13; mature T-cell proliferation-1 type B1; mature T-cell proliferation 1		

For research use only.

GENE DESCRIPTION	mature T-cell proliferation 1
GENE ID	4515
ORIGIN SPECIES	Human
ACCESSIONS	BC002600
VECTOR	pDONR223
NCBI SYMBOL	MTCP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcagaaggatccgtgccagaagcaagcctgtgagatacagaaatgtttacaagccaacagctacatggaatcaaagtgtcaggtgtcatccaagaactgcgtaagtgtgtgtcagatccaagggaagatctgtcgtctgttcaggatttgaaaaagaagaggaagaaaacct
aacacggaagtctgcatcaaagtaa