

RPS26-ribosomal protein S26 Gene

Reference: PTXBC002604

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA10; S26; 40S ribosomal protein S26; ribosomal protein S26		

For research use only.

GENE DESCRIPTION	ribosomal protein S26
GENE ID	6231
ORIGIN SPECIES	Human
ACCESSIONS	BC002604
VECTOR	pDONR223
NCBI SYMBOL	RPS26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaaagaaaaggaacaatggctgtgccccaaaggccggggccacgtgcagcctattcgtgcactaactgtgcccgatgcgtgcccaaggacaaagccattaaagaaatcgtcattcgaaacatagtgaggccgcagcagtcaggacatttctgaagcgagcgtcttga
tgcctatgtgcttccaagctgtatgtgaagctacattactgtgtgagttgtgcaattcacagcaaaagtagtcaggaatcgtatcgtgaagcccgaaggaccgaacacccccaccccgatttagacctgcgggtgctgccccagctccccaccaaagcccatgtaa