

FKBP1B-FK506 binding protein 1B, 12.6 kDa Gene

Reference: PTXBC002614

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP1B; peptidyl-prolyl cis-trans isomerase FKBP1B; FKBP12.6; FKBP1L; OTK4; PKBP1L; PPIase; 12.6 kDa FK506-binding protein; 12.6 kDa FKBP; FK506 binding protein 1B, 12.6 kDa; FK506-binding protein 12.6; FKBP-12.6; FKBP-1B; calstabin 2; h-FKBP-12; i		

For research use only.

GENE DESCRIPTION	FK506 binding protein 1B, 12.6 kDa
GENE ID	2281
ORIGIN SPECIES	Human
ACCESSIONS	BC002614
VECTOR	pDONR223
NCBI SYMBOL	FKBP1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtggagatcgagaccatctccccggagacgggaaggacattcccaagaaggccaaacgtgtgtgtgtgactacacaggaatgctccaaatgggaagaagttgattcatccagagacagaacaaaccttcaagttcagaattggcaaacaggaagtcataaaggttt
gaagaggggtgcagccagctgggtcctctttctcctctcccatctgccccatcctgctag