

SNRPE-small nuclear ribonucleoprotein polypeptide E Gene

Reference: PTXBC002639

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HYPT11; SME; Sm-E; snRNP-E; small nuclear ribonucleoprotein E; sm protein E; small nuclear ribonucleoprotein polypeptide E		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein polypeptide E
GENE ID	6635
ORIGIN SPECIES	Human
ACCESSIONS	BC002639
VECTOR	pDONR223
NCBI SYMBOL	SNRPE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtaccgtggccagggtcagaaagtcagaagggttatgggtcagcccatcaacctcatcttcagatacttacaagaatagatcgcggattcaggtgtggtctctatgagcaagtgaatatcggatagaaggctgtatcattggtttgatgatatgaacctgtattagatgatgcagaa
gagattcattctaaaacaaagtcaagaaaacaactgggtcggatcatgctaaaaggagataatattactctgtacaaagtgtctccaactag