



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

UCN2-urocortin 2 Gene

Reference: PTXBC002647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SRP; UCN-II; UCNI; URP; urocortin-2; prepro-urocortin 2; stresscopin-related peptide; ucn II; urocortin II; urocortin-related peptide; urocortin 2		

For research use only.

GENE DESCRIPTION	urocortin 2
GENE ID	90226
ORIGIN SPECIES	Human
ACCESSIONS	BC002647
VECTOR	pDONR223
NCBI SYMBOL	UCN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacacgggtgctcgtcgtgctgctgatggctgctgagctgtgggcagagctcctggtgctccagctgaccctatcccaacctccagctcgcgcctcagaattctcccgaccacacftcccgcacctgcgccctcagagagccccctcagctgctcccaatgccgctggtgggctgcccagagcca
ctgcagcccccaccgcacccctgctcgcgcattgctctatcctgctgagctcccatcgccctcttgcagatctactggagcaagcccgccggccaggctgccaggagcaggccaccacccaagcccccgcattctgcccctgtctgcccactctga