



Reference: PTXBC002649

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1.2; H1C; H1F2; H1s-1; histone H1.2; H1 histone family, member 2; histone 1, H1c; histone H1c; histone H1d; histone H1s-1; histone cluster 1, H1c; histone cluster 1 H1 family member c		

For research use only.

GENE DESCRIPTION	histone cluster 1, H1c
GENE ID	3006
ORIGIN SPECIES	Human
ACCESSIONS	BC002649
VECTOR	pDONR223
NCBI SYMBOL	HIST1H1C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgtccgacagctgtctctccgctcccgctccgctccgctcctgcggagaaagccctctgaataagaagaagcgccgcaaaaagcgctgggggtacgcctctgaaggcctgtctgtcccccgtctcagagctcatcaccaaggctgtggccgctctaaagcgtgacgggatttctctgt
ctgctctgaaaaaagcgttggtctgccgcgcctctgatgtggagaaaaacaacgcctgatacaacttggtctcaagagccttggtgagcaagggcactctgtgtgcaaacgaaaggcaccggtgctctgtgctctttaactcaacaagaaggcagcctccgggggaagccgaaccccaag
gttaaaaagctgcggcgccgagcaaaccaactagaagaagccaggttggtggcgacccaagaagccccaagaagcgcgctgcggcgcaactccccaagaagaagcgtctgaagaacacacccaagaagaagcggcgctgctatcaacaaaagtgctgaagcccaagaatag
ccaaaggttcgcaggcccaagaagaagctgcgcaaaagctctctgaagctctggaagcccaagccgctgagcccaagctgctctcaagctctaaagaagcgcgtccgaagaagcgcgtccgaagaagcgcgtccgctatcaacaaaagtgctgaagcccaagaatag