

RPS27-ribosomal protein S27 Gene

Reference: PTXBC002658

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPS-1; MPS1; S27; 40S ribosomal protein S27; metallopan-stimulin 1; metallopanstimulin 1; ribosomal protein S27		

For research use only.

GENE DESCRIPTION	ribosomal protein S27
GENE ID	6232
ORIGIN SPECIES	Human
ACCESSIONS	BC002658
VECTOR	pDONR223
NCBI SYMBOL	RPS27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctctcgcaaggatctccttcctcctcagaagaggagaagaggaacacaaagaagcgcctggtgcagagcccaattcctactcatggatgtgaaatgccaggatgctataaaatcaccacggctcttagccatgcacaaacggtagttttgtgtgttgctgtccactg
tcctctgccagctacaggagagaaaagcaaggccttacagaaggatgttccttcaggagggaagcagcactaa