

CCDC56-coiled-coil domain containing 56 Gene

Reference: PTXBC002698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC56; COX25; HSPC009; MITRAC12; cytochrome c oxidase assembly factor 3 homolog, mitochondrial; coiled-coil domain-containing protein 56; cytochrome c oxidase assembly protein 3 homolog, mitochondrial; mitochondrial translation regulation assembly interm		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 56
GENE ID	28958
ORIGIN SPECIES	Human
ACCESSIONS	BC002698
VECTOR	pDONR223
NCBI SYMBOL	CCDC56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcttcgggagctggtgacctctggaattctaagcgtggagaggccccgttcgctcagcgtatcgaccgactcgggagaagctgacacccgagcaactgcattccatgcggcagcgaggcttgcctcagtgagcagaaggctctaccacggcgcggaacccggaacatcgtg
accggcctaggcatcggggccctggtgttggtatattatggttacaccttctactcgaattccagaggcgtttcctagatgagctagaagacgaggccaaagctgcccagagccgagctctgccaaggcgtcagggctcctaa